



TECH BULLETIN

039
03/85

Published regarding engineering changes and improvements

SUBJECT: Flat Valves for Model 390

Beginning with serial number 2850101, the flat valve assembly became standard in the Model 390 pumps.

The Flat Valve parts are not interchangeable with the Quiet Valve parts and must be replaced as a set.

	OPTIONAL	STANDARD
	[Quiet Valve]	[Flat Valve]
Valve.....	43721	104302
Seat.....	43779	103956
Retainer	43780	103957
Spring	43751	43751
O-Ring, D.V.S.	43793	43793
Valve Kit.....	30859	34010

CAT PUMPS
Technical Services Department

03/96